

VetScan
VS2 Chemistry Analyzer



ABAXIS

BETTER + BETTER



On-Site Diagnostic Testing, with 10 Panels to Choose From

The VetScan VS2 is a state-of-the-art chemistry, electrolyte, immunoassay and blood gas analyzer that delivers uncompromising accuracy from just two drops of whole blood, serum or plasma. The VetScan VS2 is amazingly simple and intuitive, featuring an advanced user interface, expanded printing capabilities and precision reference laboratory quality results in 12 minutes—making it ideal for veterinary clinics, mobile practitioners, research laboratories, and pharmaceutical and biotech companies.

Features

- Chemistry, electrolyte, immunoassay and blood gas with one simple instrument

- Accuracy comparable to reference laboratories

- Speed: 12 minutes to results

- Intelligent Quality Control (iQC)

- Sample integrity alerts

- 100 μ L sample size – whole blood, serum or plasma

- Cost effective

- Portable, space-saving design

- Set and archive reference ranges

- Stores 5,000 patient records

- Connectivity

Benefits

- Save time, counter space and money while improving in-house testing with the availability of T₄, Bile Acids, tCO₂ and Canine Heartworm

- Offers peace of mind to diagnose and treat patients appropriately

- Obtain results, diagnose, discuss with clients and begin treatment within one visit

- iQC monitors the analyzer, rotor and sample before and during every analysis to ensure reliable results

- Monitors interference to ensure accurate results and reports levels of HEM, LIP and ICT to provide useful information on patients and sample handling

- Test a full chemistry profile using only two drops of sample from multiple species – vital in smaller patients


- Reduce your costs and save time:
 - Maintenance-free
 - No sample preparation – test whole blood

- The 11-pound analyzer has the ability to fit in one square foot of space and can be used anywhere and anytime

- Enter pre-labeled or unique species reference intervals as well as archive ranges for later use

- Provides peace of mind that results are secure and allows recall of patient history

- Connect with Practice Management Software, other VetScan analyzers, external printers or a keyboard using the four USB ports



Comprehensive Test Portfolio – 10 Different Profiles to Choose From

Analytes	ALB	ALP	ALT	AMY	AST	BA	BUN	Ca	CHOL	CHW	CK	CL-	CRE	GGT	GLOB*	GLU	K+	Mg	Na+	PHOS	T ₄	TBIL	tCO ₂	TP	UA
Canine Wellness Profile† #500-0044	•	•	•				•	•		•			•		•	•				•		•		•	
Comprehensive Diagnostic Profile #500-0038	•	•	•	•			•	•					•		•	•	•		•	•		•		•	
Prep Profile II #500-0026		•	•				•						•			•								•	
Critical Care Plus #500-0042			•				•					•	•			•	•		•				•		
T ₄ /Cholesterol Profile #500-0037									•												•				
Kidney Profile Plus† #500-0045	•						•	•				•	•			•	•		•	•			•		
Mammalian Liver Profile #500-0040	•	•	•			•	•		•					•									•		
Equine Profile Plus #500-0043	•				•		•	•			•		•	•	•	•	•		•			•	•	•	
Avian/Reptilian Profile Plus #500-0041	•				•	•	•				•				•	•	•		•	•				•	•
Large Animal Profile #500-0023	•	•			•		•	•			•			•	•			•		•				•	

*Calculated value.

†For VetScan VS2 only.

Intelligent Quality Control (iQC)

With the VetScan VS2 Chemistry Analyzer quality control is built-in—thanks to a sophisticated iQC system that combines with proprietary algorithms to assure quality dependability with every result.



Simple Operation

Practices everywhere have discovered how the VetScan Chemistry Analyzer has improved practice efficiencies and saved them time and money. Simplicity, speed, through-put with results on-site, in minutes.



1 Insert sample.



2 Insert rotor.



3 Read results.

Requires only two drops of whole blood.

VetScan VS2

Easy as 1-2-3

For more information on the VetScan VS2, contact your Abaxis Representative at 800.822.2947, email us at vetscan@abaxis.com or click www.abaxis.com/veterinary.

Abaxis is your total on-site laboratory partner, offering the comprehensive and innovative VetScan line of diagnostic instruments and individual tests. The VetScan name has become synonymous with convenience, ease of use and cost-effectiveness, delivering on our promise of speed, accuracy and reliability. It's not just better technology. It's a better way.



**Abaxis, Inc.
Worldwide Headquarters**

3240 Whipple Road
Union City, CA 94587
Tel 800.822.2947
Fax 510.441.6150
Email vetscan@abaxis.com

Abaxis Europe GmbH

Otto-Hesse-Strasse 19
T9, 3. OG Ost
D-64293 Darmstadt Germany
Tel +49 6151 350 790
Fax +49 6151 3507911
Email abaxis@abaxis.de

www.abaxis.com