

### Medical Suction Unit



**\$160**  
**\$250 with Rechargeable Battery**

Suction Units are designed for portable operation in the effective use of aspirating fluids from the airways of patients by trained, qualified medical personnel. This suction unit is ideal for providing effective aspiration in both the institutional and homecare environments. Choose from two convenient Models, the MIG-8100DC Model or the MIG-8110 A/C, D/C Model that includes a rechargeable Lithium LiFePO4 Battery.

- Effective Removal of Respiratory Fluids
- Removes Material from Wounds
- Durable, Portable System
- Light Weight DC or AC Operation
- Replaceable In-Line Bacteria Filter
- Removal Fill Tank
- Two Year Warranty

### Compact Compressor Nebulizer System



**\$70 With pediatric and adult mask**  
**Fire & Rescue Nebulizer System**

Contains Standard Circuit & NebJet  
Reusable High Capacity JetNeb  
Attractive Sleek Design  
Durable, Quiet System  
Piston Driven Compressor  
Integrated Handle for Transport  
Exclusive Lifetime Warranty



**\$55**

Child friendly design  
Lifetime Warranty!  
Weighs only 4.3 lbs  
Medication capacity 5ml  
Your active little one will enjoy the eye-catching and interactive Fire & Rescue Nebulizer System by. The Fire & Rescue Nebulizer System was designed to engage the user and ensure the proper treatment is administered. Weight: 4.3 lbs

Fire & Rescue Nebulizer System contents:  
Pediatric Compressor  
Nebulizer Cup  
Tubing  
Pediatric Mask  
Mouthpiece  
5 filters

### Interactive Building Block Compressor Nebulizer



**\$55 , Kit \$11**

Child friendly  
Eliminates anxiety, increases coordination  
Weighs only 4.3 pounds  
Lifetime Warranty  
This pediatric nebulizer engages the user's imagination while increasing the hand and eye coordination. The child will be at ease and consumed while receiving the required treatment. Parents will appreciate the interactive options and convenience the building Block feature will offer as well as the technology.

Interactive Building Block Compressor Nebulizer contents :  
Pediatric Compressor  
Nebulizer Cup  
Tubing  
Pediatric Mask  
5 Filters  
Carrying case

Options: Kit Heliport, Kit Pirates Island

### Panda Pink Compressor System



**\$50**

Esthetically appealing  
Simple to use  
Pink carrying case  
Lifetime Warranty  
The Panda Compressor System brings appeal and ease of use together making treatments child friendly and uncomplicated. The esthetic appeal of the panda compressor invites the child in for closer observation. Caregivers will find treatments are easier to administer when the child is comfortable with the compressor. The Panda Compressor is easy to transport in the carrying case provided that also holds all the necessary accessories for treatments.

Panda Pink Pediatric Compressor Nebulizer System contents :  
Pediatric Compressor  
Nebulizer Cup  
Tubing  
Pediatric Mask  
5 Filters  
Carrying case

### Panda Compressor System



**\$50**

Child friendly  
Eliminates anxiety, increases coordination  
Simple to use  
Lifetime Warranty

Panda Pediatric Compressor Nebulizer System contents :  
Pediatric Compressor  
Nebulizer Cup  
Tubing  
Pediatric Mask  
5 Filters  
Carrying case

### **Penguin Compressor Nebulizer**



**\$50**

Child friendly  
Eliminates anxiety, increases coordination  
Simple to use  
Lifetime Warranty

Penguin Compressor Nebulizer System contents :  
Pediatric Compressor  
Nebulizer Cup  
Tubing  
Pediatric Mask  
5 Filters  
Carrying case

### **Elephant - Puppy Ultrasonic Cool Mist Humidifier**



**\$50 each**

Durable Ultrasonic System  
Whisper Quiet Technology  
Variable Humidity Control  
Auto-Off Function  
Easy Fill Tank  
One Year Warranty

## View SpO2 Finger Tip Pulse Oximeter



**\$69**

### Features:

Dual color HIGH RESOLUTION  
Easy to Read OLED display.  
Display SpO2, Pulse Rate ,  
Features Pulse Bar and Plethysmograph.  
Multidirectional (Six) display modes for easy reading from all sides.  
Automatic power off with battery-low indicator.  
Normal battery operation is approximately thirty hours.  
Supports 2 AAA alkaline or rechargeable batteries.  
Ten adjustable brightness settings.  
Accommodates wide range of finger sizes from pediatric to adult.  
Sturdy lightweight design ( total 0.06 lbs without batteries ).  
Two year warranty.

### Specifications

#### SpO2

Measurement range 70-100% Resolution 1%  
Measurement accuracy 80%-100%;  $\pm 2\%$  70%-79% ;  $\pm 3\%$  0%-69 %;unspecified

#### Pulse Rate

Measurement range 30-235 bpm Resolution 1bpm  
Measurement accuracy 30 - 100;  $\pm 2$  bpm 101- 235;  $\pm 2\%$

#### Environmental

Operation Temperature 41°F to 104°F Storage Temperature -4°F to 131°F Operation Humidity RH -  $\leq 85\%$  Storage Humidity RH -  $\leq 85\%$

#### Electromagnetic Compatibility

This product complies with IEC60601-1-2 for electromagnetic compatibility Class B

#### Display

*Type:* Dual color High Resolution OLED display  
*Parameters:* SpO2, Pulse Rate, Pulse Bar and Plethysmogram  
*Mode:* 6 display modes for easy reading from all sides  
*Brightness:* Adjustable 10 levels  
Interference Resistance Capacity against Ambient Light  
Deviation is smaller than  $\pm 1\%$  between values of Oxyhemoglobin measured in natural lighting indoor condition and present lighting sources and that measured in dark room.  
Mechanical  
Dimension 2.3 in(W) x 1.3 in(H) x 1.3 in(D) Weight 0.06 lbs(without batteries)

## Fingertip Pulse Oximeter



**\$45**

The MIG-8300 is intended for noninvasive spot-check oxygen saturation measurements of the arterial hemoglobin (SpO2). This product is ideal for both adult and pediatric patients.

### Product Features Include

- SpO2%, and Pulse Rate
- Large LED Display
- Easy to Use
- Tests both Pediatric and Adult
- Light, Compact Design
- Pocket Size Convenience
- Weight is 50g including batteries.
- Powered by two AAA batteries
- Low Battery indicator
- Auto Power Off